

9853 Pacific Heights Blvd. Suite D. San Diego, CA 92121, USA Tel: 858-263-4982

Email: info@abeomics.com

## 32-13519: ProSpec 7-Band Marker

Alternative Name: Protein Marker, Protein Ladder, Protein Standard, Mw Marker, Mw ladder, Mw Standard

## **Description**

Source: E.Coli

Sterile Filtered blue solution.

ProSpec 7-Protein marker containing 8.5, 13.5, 18, 28, 40, 57 & 70 kDa. The recombinant proteins are produced in E.Coli is and purified by proprietary chromatographic techniques.

## **Product Info**

**Amount:** 0.25 ml / 0.5 ml

Content: ProSpec 7-protein marker solution (0.5ug/ul/protein) contains 12mM Tris-HCl pH-6.8, 5%

glycerol, 0.4% SDS 2.88mM 2-mercaptoethanol & 0.02% bromophenol blue.

Storage condition:

Storage condition:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Available traces that guides.

of time. Avoid multiple freeze-thaw cycles.